

KNOCHENBOHRER

BONE TWIST DRILLS

AO-ANSCHLUSS -DREISCHNEIDIG AO-SHANKS -THREE-EDGED	Durchmesser in mm	Gesamtlänge Total Length in mm	Arbeitslänge in mm	Spirallänge in mm
KG.049.25	Ø 2,5	230	205	30
KG.049.32	Ø 3,2	230	205	30
KG.049.45	Ø 4,5	195	170	55



AO-ANSCHLUSS AO-SHANKS				
KG.050.01	Ø 1,1	60	35	20
KG.050.10	Ø 1,5	85	60	23
KG.050.12	Ø 1,5	110	85	50
KG.050.20	Ø 2,0	100	75	50
KG.050.22	Ø 2,0	125	100	50
KG.050.30	Ø 2,5	105	80	50
KG.050.32	Ø 2,5	125	100	50
KG.050.36	Ø 2,5	180	155	50
KG.050.40	Ø 2,7	100	75	50
KG.050.42	Ø 2,7	125	100	50
KG.050.53	Ø 3,2	145	120	50
KG.050.57	Ø 3,2	195	170	50
KG.050.62	Ø 3,5	110	85	50
KG.050.67	Ø 3,5	195	170	50
KG.050.73	Ø 4,5	145	120	50
KG.050.77	Ø 4,5	195	170	50
KG.050.87	Ø 6,0	195	170	60



BOHRER RUNDSCHAFT DRILLS ROUND SHANK	Durchmesser in mm	Gesamtlänge Total Length in mm
KG.055.10	Ø 1,0	130
KG.055.15	Ø 1,5	130
KG.055.20	Ø 2,0	130
KG.055.25	Ø 2,5	130
KG.055.27	Ø 2,7	130
KG.055.32	Ø 3,2	130
KG.055.35	Ø 3,5	130
KG.055.45	Ø 4,5	130
KG.055.60	Ø 6,0	130

